

pluriSpin® Human CD19⁺ Cell Depletion (10 ml)

Order No.:19-06002-32



Description

pluriSpin® Human CD19⁺ Cell Depletion is designed to isolate pure, viable and untouched b cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels B cells, while the isolated human b cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified b cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|----------------------|-------------------------------------------------------------------------------|
| Size | 10 ml |
| Species | Human |
| Cell Types | B cells |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

