

1x PBS Solution (10 * 500 ml)

Order No.:60-00000-12



Description

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

Components

- NaCl (136,89 mMol/L)
- Na₂ HPO₄ • 2H₂ O (10mMol/L)
- KH₂ PO₄ (1,98 mMol/L)
- KCl (2,68 mMol/L)

Additional Information

Size	10 * 500 ml
Delivery Time (days)	1-2
Area of interest	Cell Biology
Sterility	Sterile
pH Range	pH 7.4
Shipping Condition	Room Temperature

Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.