

universal M-pluriBead® anti-rat

Order No.:31-GTaRT-20



Description

pluriBead universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

Additional Information

| | |
|----------------------|-------------------------------------------------------------------------------|
| Size | 2 ml |
| Species | Rat (Rattus) |
| Cell Types | Any |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.