

CD55 M-pluriBead® anti-hu

Order No.:19-05500-20



Description

CD55 pluriBead® suspension for the positive cell isolation of CD55+ target cells from any human sample material. (2ml)

Additional Information

| | |
|----------------------|---|
| Size | 2 ml |
| Species | Human |
| Cell Types | Hematopoietic cells, Non-hematopoietic cells |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.