

CD271 M-pluriBead® anti-hu

Order No.:19-27100-20



Description

CD271 pluriBead® suspension for the positive cell isolation of CD271+ target cells from any human sample material. (2ml)

Additional Information

| | |
|----------------------|---|
| Size | 2 ml |
| Species | Human |
| Cell Types | Mesemchymal cells, Schwann cells, Stem/Stromal cells |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.