

F4/80 S-pluriBead® anti-ms

Order No.:29-F4800-10



Description

F4/80 pluriBead® suspension for the positive cell isolation of F4/80+ target cells from any murine sample material. (2ml)

Additional Information

Size	2 ml
Species	Mouse (Murine)
Cell Types	Monocytes, Macrophages
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.