

## Funnel (50 pcs. , non-sterile)



**Order No.:**42-50000-50

### Description

Funnel - Sterile device to support sample load onto pluriStrainer.

### Additional Information

Size	50 pcs.
Delivery Time (days)	1-2
Variation	50 pcs, non-sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.