

pluriMag Rack I (1.5 / 2 ml tubes)

Order No.:50-01010-01



Description

pluriMag Rack is a magnetic holder for small-scale protein and nucleic acids purification and cell enrichment or depletion with magnetic beads. pluriMag Rack can be used for sample volumes of up to 5 ml. The holder make an excellent complement to the following pluriMag products (containing magnetic beads) from pluriSelect. pluriMag Rack works also with all standard magnetic micro beads.

It is optimized for efficient magnetic separation of all types of magnetic (micro) beads (range 1–8 μm in diameter) from the supernatant (10–2,000 μl) in seconds, allowing for easy pipette removal without sample loss. The pluriMag Rack combines a strong magnetic field for efficient separations in standard microcentrifuge tubes.

- Best working volume: 10–2,000 μl (1ml -50ml)
- Holds up to 10 standard 1.5–2 ml microcentrifuge tubes in or 4 flow cytometry tubes
- Good control and visibility of your sample

Additional Information

| | |
|----------------------|---|
| Size | 1 device |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.