

Flow Cytometry Tube 5 ml, sterile



	Order No.	Packaging	Variant
<input type="checkbox"/>	05-03040-70	20 Bags, 25 pcs / bag	500 pcs, sterile
<input type="checkbox"/>	05-03040-40	1 Bag, 25 pcs / Bag	25 pcs, sterile
<input type="checkbox"/>	05-03040-60	4 Bags, 25 pcs / bag	100 pcs, sterile
<input type="checkbox"/>	05-03040-50	2 Bags, 25 pcs / bag	50 pcs, sterile

Description

5 ml round-bottom tubes are suitable for a wide range of laboratory applications such as sample preparation or sample storage. The tubes are ideal for flow cytometry. The bundle of pluriStrainer Mini and flow cytometry tubes can be used for sample preparation and as replacement for Falcon® round-bottom tubes with blue cell strainer cap (35 µm nylon mesh). Use the pluriStrainer Mini prior fluorescence-activated cell sorting (FACS) to get a single cell suspension and avoid clogging the flow cell.

Application

- Sample preparation, storage and centrifugation
- Cell isolation and enrichment
- Cell staining for flow cytometry analysis

Additional Information

Prefilled	Unfilled
Delivery Time (days)	1-2
Storage Condition	Room Temperature
Shipping Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

