

CD11b S-pluriBead® anti-ms

Order No.:29-011b0-10



Description

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any murine sample material. (2ml)

Additional Information

| | |
|----------------------|---|
| Size | 2 ml |
| Species | Mouse (Murine) |
| Cell Types | Monocytes, Granulocytes, Macrophages, NK cells (subset) |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.