

pluriSpin® Human Monocyte Enrichment (10 ml)

Order No.:19-01001-32



Description

pluriSpin® Human Monocyte Enrichment is designed to isolate pure, viable and untouched monocytes from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, B cells, NK cells and granulocytes, while the negatively isolated human monocytes are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|----------------------|---|
| Size | 10 ml |
| Species | Human |
| Cell Types | Monocytes |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

