

## universal S-pluriBead® anti-ms

**Order No.:**31-GTaMS-10



### Description

pluriBead universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

### Additional Information

|                      |   |
|----------------------|---|
| Size                 | 2 ml  |
| Species              | Mouse (Murine)  |
| Cell Types           | Any   |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >90%  |
| Storage Condition    | 4-8 °C  |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.