

## pluriSpin® Human CD8<sup>+</sup> T Cell Enrichment



Order No.	Size
19-04001-31	2 ml
19-04001-32	10 ml

### Description

pluriSpin® Human CD8<sup>+</sup> T Cell Enrichment is designed to isolate pure, viable and untouched cytotoxic T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, helper T cells, B cells, NK cells and granulocytes, while the negatively isolated human cytotoxic T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic T cells are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Species	Human
Cell Types	Cytotoxic T cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.