


Centrifuge Tube 50 ml (sterile, 300 pcs.)



Order No.: 05-00001-01

Description

50 ml tubes are suitable for a wide range of laboratory applications such as sample preparation or sample storage. The tubes are ideal for the use with our **pluriStrainers** and guarantee a perfect fit with the Connector Ring. Moreover the tubes can be used as container for **density gradient centrifugation** or the enrichment of primary cells with **pluriBead**.

Screw tube, 50 ml, 114 x 28 mm, pointed bottom, PP, with printed writing field and scaling, with mounted red cap, sterile and pyrogen-free / endotoxin-free, 25 pcs./bag

Application

- Mixing
- Incubation
- Centrifugation
- Storage
- Collection vessel for filtration and purification of liquids

Additional Information

Size	300 pcs.
Prefilled	Unfilled
Delivery Time (days)	1-2
Variation	20 pcs, sterile
Packaging	25 pcs. / Bag
Storage Condition	Room Temperature
Shipping Condition	Room Temperature

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.