

pluriSpin® Human CD4⁺ T Cell Enrichment



| Order No. | Size |
|-------------|-------|
| 19-03001-31 | 2 ml |
| 19-03001-32 | 10 ml |

Description

pluriSpin® Human CD4⁺ T Cell Enrichment is designed to isolate pure, viable and untouched CD4⁺ T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, cytotoxic T cells, B cells, NK cells and granulocytes, while the negatively isolated human CD4⁺ T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD4⁺ T cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|----------------------|---|
| Species | Human |
| Cell Types | Helper T cells |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.