

## ÜberStrainer 1 µm



**Order No.:**43-70001-03

### Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

### Additional Information

|                      |                                                               |
|----------------------|---------------------------------------------------------------|
| Mesh Size            | 1 µm                                                          |
| Size                 | 10 pcs                                                        |
| Fabric material      | PET (Polyethylenterephthalat)                                 |
| Housing Material     | PP (Polypropylene)                                            |
| Mesh Type            | woven                                                         |
| Sterility            | Sterile                                                       |
| Storage Condition    | Room Temperature                                              |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

