

ÜberStrainer 400 µm



Order No.:43-70400-03

Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

Additional Information

| | |
|----------------------|---|
| Mesh Size | 400 µm |
| Size | 10 pcs |
| Fabric material | PET (Polyethylenterephthalat) |
| Housing Material | PP (Polypropylene) |
| Mesh Type | woven |
| Sterility | Sterile |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

