

## pluriSpin® Human PBMC Depletion (10 ml)

Order No.:19-08002-32



### Description

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

### Additional Information

|                      |   |
|----------------------|---|
| Size                 | 10 ml   |
| Species              | Human   |
| Cell Types           | Non-PBMC cells  |
| Storage Condition    | 4-8 °C  |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.