

## Monocytes Spin Medium, 100 ml

**Order No.:**60-00095-10



### Description

Monocytes Spin Medium is a ready-to-use, sterile density medium for the enrichment of monocytes in a high yield from fresh peripheral blood or buffy coat through single step density centrifugation. Depending on the sample a purity of 60-80% is possible. Monocytes Spin is a gentle alternative to bead-mediated monocyte enrichments, especially for downstream applications like cell culture and differentiation.

### Additional Information

|                      |  |
|----------------------|--|
| Size                 | 100 ml   |
| Delivery Time (days) | 1-2  |
| Storage Condition    | Room temperature, protect from light                         |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range             | 6.8 - 7.4  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.