

## 10x PBS Solution (500 ml)

**Order No.:**60-00010-11



### Description

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- RNase/DNase free
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (10mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

Size	500 ml
Sterility	Sterile
pH Range	pH 6.8
Shipping Condition	Room Temperature
Storage Condition	Room Temperature

**PluriSelect USA**

2030 Gillespie Way, Unit 103, CA 92020 El Cajon, US  
 Phone: +1 619-202-4297 | Fax: +1 619-828-8003  
 sales-support@pluriselect-usa.com  
<https://www.pluriselect.com/>

Regulatory Statement

For research use only. Not for use in diagnostic procedures,

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.