

## 1x PBS Solution (500 ml)

**Order No.:** 60-00000-10



### Description

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without**:

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (10mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

Size	500 ml
Delivery Time (days)	1-2
Area of interest	Cell Biology
Sterility	Sterile
pH Range	pH 7.4
Shipping Condition	Room Temperature

#### PluriSelect USA

2030 Gillespie Way, Unit 103, CA 92020 El Cajon, US  
 Phone: +1 619-202-4297 | Fax: +1 619-828-8003  
[sales-support@pluriselect-usa.com](mailto:sales-support@pluriselect-usa.com)  
<https://www.pluriselect.com/>

Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.