

## CD31 M-pluriBead® anti-ms



**Order No.:**29-03100-20

### Description

CD31 pluriBead® suspension for the positive cell isolation of CD31+ target cells from any murine sample material. (2ml)

### Additional Information

Size	2 ml
Species	Mouse (Murine)
Cell Types	Platelets, Endothelial cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.