

## CD14 M-pluriBead® anti-porcine

**Order No.:**23-01400-20



### Description

CD14 pluriBead® suspension for the positive cell isolation of CD14 target cells from any porcine sample material. (2ml)

### Additional Information

Size	2 ml
Species	Pig (Porcine)
Cell Types	Monocytes
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.