

CD62L M-pluriBead® anti-porcine

Order No.:23-062L0-20



Description

CD62L pluriBead® suspension for the positive cell isolation of CD62L+ target cells from any porcine sample material. (2ml)

Additional Information

Size	2 ml
Species	Pig (Porcine)
Cell Types	Monocytes, T cells, Granulocytes, B cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.