

CD8b M-pluriBead® anti-rat

Order No.:22-008b0-20



Description

CD8b pluriBead® suspension for the positive cell isolation of CD8b+ target cells from any rat sample material. (2ml)

Additional Information

Size	2 ml
Species	Rat (Rattus)
Cell Types	Cytotoxic T cells, Thymocytes
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.