

## CD271 S-pluriBead® anti-sheep

**Order No.:**24-27100-10



### Description

CD271 pluriBead® suspension for the positive cell isolation of CD271+ target cells from any sheep sample material. (2ml)

### Additional Information

Size	2 ml
Species	Sheep
Cell Types	Mesemchymal cells, Schwann cells, Stem/Stromal cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.