

## CD9 S-pluriBead® anti-bovine

**Order No.:**25-00900-10



### Description

CD9 pluriBead® suspension for the positive cell isolation of CD9+ target cells from any bovine sample material. (2ml)

### Additional Information

Size	2 ml
Species	Cow (Bovine)
Cell Types	Platelets, Granulocytes
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.