

10% BSA Stock Solution in 1x PBS


☐
☐
☐

Order No.	Size	Sterility
60-00020-10	50 ml	Sterile
60-00020-11	250 ml	Sterile
60-00020-12	500 ml	Sterile

Description

10% BSA stock solution in 1x PBS can be used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

Additional Information

Delivery Time (days)	1-2
Area of interest	Cell Biology
Sterility	Sterile
pH Range	pH 7.4
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

