

**Funnel (50 pcs. , non-sterile)****Order No.:**42-50000-50**Description**

Funnel - Sterile device to support sample load onto pluriStrainer.

**Additional Information**

Size	50 pcs.
Delivery Time (days)	1-2
Variation	50 pcs, non-sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

**Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.