

Funnel


☐
☐
Order No.
Variant

42-50000-03

25 pcs, sterile

42-50000-50

50 pcs, non-sterile

Description

Funnel - Sterile device to support sample load onto pluriStrainer.

Additional Information

Size	50 pcs.
Delivery Time (days)	1-2
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.