

## Funnel



**Order No.**
**Variant**

42-50000-03

25 pcs., sterile

42-50000-50

50 pcs., non-sterile

### Description

Funnel - Sterile device to support sample load onto pluriStrainer.

### Additional Information

Size	50 pcs.
Delivery Time (days)	1-2
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.