

## Funnel


☐
☐
**Order No.**
**Variant**

42-50000-03

25 pcs, sterile

42-50000-50

50 pcs, non-sterile

### Description

Funnel - Sterile device to support sample load onto pluriStrainer.

### Additional Information

|                      |  |
|----------------------|--|
| Size                 | 50 pcs.  |
| Delivery Time (days) | 1-2  |
| Shipping Condition   | Room Temperature   |
| Storage Condition    | Room Temperature   |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.