

Syringe Strainer 100 µm (12 pcs., non-sterile)

Order No.:43-71100-50



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

Size	12 pcs.
Flow Control	Yes
Mesh Size	100 µm
Housing Material	LD-PE (Low Density Polyethylen)
Mesh Type	woven
Sterility	Non-Sterile
Variation	12 pcs, non-sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.