

## Syringe Strainer 20 µm


☐
☐

Order No.	Sterility	Variant
43-71020-50	Non-Sterile	12 pcs, non-sterile
43-71020-51	Sterile	12 pcs sterile

### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Mesh Size	20 µm
Size	12 pcs.
Flow Control	Yes
Housing Material	LD-PE (Low Density Polyethylen)
Mesh Type	woven
Sterility	Sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

