

Syringe Strainer 30 µm


☐
☐

Order No.	Sterility	Variant
43-71030-50	Non-Sterile	12 pcs, non-sterile
43-71030-51	Sterile	12 pcs sterile

Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

Mesh Size	30 µm
Size	12 pcs.
Flow Control	Yes
Housing Material	LD-PE (Low Density Polyethylen)
Mesh Type	woven
Sterility	Sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

