

Syringe Strainer 50 µm


☐
☐

| Order No. | Sterility | Variant |
|-------------|-------------|---------------------|
| 43-71050-50 | Non-Sterile | 12 pcs, non-sterile |
| 43-71050-51 | Sterile | 12 pcs sterile |

Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

| | |
|----------------------|--|
| Mesh Size | 50 µm |
| Size | 12 pcs. |
| Flow Control | Yes |
| Housing Material | LD-PE (Low Density Polyethylen) |
| Mesh Type | woven |
| Sterility | Sterile |
| Shipping Condition | Room Temperature |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

