

## Syringe Strainer 85 µm


☐
☐

| Order No.   | Sterility   | Variant             |
|-------------|-------------|---------------------|
| 43-71085-50 | Non-Sterile | 12 pcs, non-sterile |
| 43-71085-51 | Sterile     | 12 pcs sterile      |

### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

|                      |  |
|----------------------|--|
| Mesh Size            | 85 µm  |
| Size                 | 12 pcs.  |
| Flow Control         | Yes  |
| Housing Material     | LD-PE (Low Density Polyethylen)                              |
| Mesh Type            | woven  |
| Sterility            | Sterile  |
| Shipping Condition   | Room Temperature   |
| Storage Condition    | Room Temperature   |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

