

pluriSpin® Human CD4⁺ T Cell Enrichment (2 ml)

Order No.:19-03001-31



Description

pluriSpin® Human CD4⁺ T Cell Enrichment is designed to isolate pure, viable and untouched CD4⁺ T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, cytotoxic T cells, B cells, NK cells and granulocytes, while the negatively isolated human CD4⁺ T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD4⁺ T cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Size	2 ml
Species	Human
Cell Types	Helper T cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

