

pluriSpin® Human CD4⁺ Cell Depletion (10 ml)

Order No.:19-03002-32



Description

pluriSpin® Human CD4⁺ Cell Depletion is designed to isolate pure, viable and untouched helper t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Helper T cells, while the isolated human helper t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified helper t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Size	10 ml
Species	Human
Cell Types	Helper T cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

