

## pluriSpin® Human T Cell Depletion





Order No.	Size
19-02002-31	2 ml
19-02002-32	10 ml

### Description

pluriSpin® Human T Cell Depletion is designed to isolate pure, viable and untouched t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, while the isolated human t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Delivery Time (days)	1-2
Species	Human
Cell Types	T cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

